

Lifblood Shippers - Receipt and Use by External Institutions

Purpose

This Work Instruction (WI) describes the process to:

- Provide information on the shipper packing configuration used by Lifeblood to transport blood products and actions you should take on receipt of a shipper from Lifeblood.
- Provide instructions to non-Lifeblood institutions that use Lifeblood shippers for their own purpose either internally or externally.

Disclaimer

This document shows the packing configurations used by Lifeblood staff when packing Lifblood shippers for transporting blood products to customers (Approved Health Providers). While Lifeblood considers that these procedures are appropriate for its own purposes, Lifeblood does not acknowledge, or have control of the particular circumstances of Approved Health Providers

Your attention is drawn to the following;

- *Customers must consider the transport requirements of Blood Products, the transport mechanism and environmental aspects in determining the appropriate packing configurations.*
- *Customers should note that Lifeblood strongly recommends the use of temperature recording data loggers for transportation over long distances and/or where extended transit times may be experienced, as the shippers are only validated for the transport of blood products over distances or timeframes according to Lifeblood specifications contained herein.*

Lifeblood cannot accept responsibility for the viability of blood products transported in a Lifblood shipper that was not packed by Lifeblood staff.

Parent procedure

This WI forms part of the procedure described in *Issuing and Packing components and Products* (SOP-00326).

Before you begin

N/A

1. Receipt of shippers

The Lifblood National Shipper Configurations ensure that blood and blood products remain within the required temperature specification during transportation. There are several packing options, the use of which will depend on:

- The component type
- The ambient temperature
- The number of components
- The anticipated transit time

Please ensure all shipments received from Lifeblood meet the criteria stated below.

If...	then...
the shipper or bag appears to have been tampered	<ul style="list-style-type: none"> • quarantine the units • contact your nearest Lifeblood Customer Service Department as soon as possible.

Lifeblood Shippers - Receipt and Use by External Institutions

If...	then...
the shipper configuration does not meet the packing criteria	<ul style="list-style-type: none"> take a manual surface temperature of the contents by placing a thermometer between two units and leave for 2-5 minutes to allow the temperature to stabilise quarantine the units contact your nearest Lifeblood Customer Service Department as soon as possible.

2. Preparation of shippers

Always check that the shippers are clean, inside and out, before use. Where there are any visible signs of damage or contamination with blood residue, the packing material is to be discarded. Deface and/or remove all old labels on the outer carton prior to discard.

2.1 Material required for shipper use

- Single insulated shipper (with foam inserts)
- Room temperature ballast packs (conditioned at +20°C to +24°C for 24hrs prior to use)
- Chilled ballast packs (conditioned at +2°C to +6°C for 24hrs prior to use)
- Frozen ballast packs (conditioned at approximately -19°C for 48hrs prior to use)
- Plastic liner bag
- Foil pouch
- Cardboard dividers with holes
- Data Loggers (data loggers for customer’s own use to be supplied by customer)

NOTE: Please refer to Appendix A for supplier details

2.2 Storage of ballast material

Table 1. Pre-conditioning and storage times for ballast packs.

Store ballast packs under the relevant conditions for the minimum time as described below.			
Ballast type	Pre-conditioned temperature	Minimum time	Configurations
Room temperature	+20°C to +24°C	24hrs	P1, P2
Chilled	Approximately +4°C	24hrs	R1, R2, R3, R4
Frozen-19°C	Approximately -19°C (domestic freezer)	48hrs	R1, R4
Frozen -40°C	Approximately -40°C	24hrs	F1
<p><u>NOTE:</u> Where the use of a data logger is necessary, the data logger should be conditioned for the temperature range required in accordance with manufacturer’s instructions.</p> <p>Failure to condition ballast packs appropriately may adversely affect internal shipper temperature.</p>			

2.3 Ambient temperature range

All configurations herein were validated using an external ambient temperature range of +4°C to +42°C.

Lifeblood Shippers - Receipt and Use by External Institutions

3. Red cells

Red cell components should be transported using one of **Packing Configurations R1, R2, R3 or R4**, the choice of which will depend on the anticipated transport time. If the transport time is estimated to approach or exceed the validated time, shown below, or if product is transported by independent courier (e.g. bus or aircraft), a data logger must be included.

Configurations R1, R2, R3 and R4 are designed to maintain the components within a temperature range of 2°C to 10°C as recommended in the Council of Europe "Guide to the Preparation, Use and Quality Assurance of Blood Components", 14th edition.

Table 2 - Packing configurations and the validated transport time.

Packing Configuration	No. of components per shipper	Validated Transport time (VTT)
R1	1 to 10 red cells	6 hours** Do not use for air transit
R2	1 to 12 red cells	3 hours 25 mins*
R3	1 to 14 red cells	8 hours 25 mins*
R4	1 to 10 red cells	16 hours 18 mins*











* If anticipated transport time exceeds the maximum VTT, a data logger must be placed in the shipper between the components

**The R1 shipper configuration cannot be used for anticipated transport times exceeding it's VTT of 6 hours, even with a data logger in use






All configurations require a specific number of units to maintain temperature. Where fewer red cell units are to be packed, the specified number of chilled ballast packs (approximately +4 °C) must be added. The correct number of ballast packs for each configuration is shown in the packing diagrams below.

The ballast packs should be placed outside the plastic liner to prevent contamination of the red cells, should the ballast packs leak.







Lifeblood Shippers - Receipt and Use by External Institutions

R1	No. Components per shipper	Validated Transport Time [#]	Legend	Name	Quantity	Note
		1 to 10 red cell units	6hr		Foil pouch	1
				Frozen -19°C	2	
				Chilled ballast	2	
				Cardboard divider	1	Used to stop contact with ballast
				Extra chilled ballast (as per table) Place ballast outside the plastic liner bag		
					Number of red cell units	Number of additional ballast packs
					4 to 10	0
					3	1
					2	2
					1	3
					1 to 3 paediatric units	4
				Plastic liner bag	1	Used to contain red cells and logger if added
				Cardboard divider	1	Used to stop contact with ballast
				Chilled ballast	2	
Notes:	# Cannot be transported by air # Cannot be used for anticipated transport times exceeding its VTT, even with a data logger in use			Frozen -19°C	2	Place on bottom











Lifeblood Shippers - Receipt and Use by External Institutions

R2	No. of Components per Shipper	Validated Transport Time [#]	Legend	Name	Quantity	Note					
	1 to 12 red cell units	3hr 25min		Chilled ballast	2						
		Extra chilled ballast (as per table) Place ballast outside the plastic liner bag	<table border="1"> <thead> <tr> <th>Number of red cell units</th> <th>Number of additional ballast packs</th> </tr> </thead> <tbody> <tr> <td>1 to 12</td> <td>0</td> </tr> <tr> <td>1 to 3 paediatric units</td> <td>4</td> </tr> </tbody> </table>	Number of red cell units	Number of additional ballast packs	1 to 12	0	1 to 3 paediatric units	4		
	Number of red cell units	Number of additional ballast packs									
	1 to 12	0									
1 to 3 paediatric units	4										
	Plastic liner bag	1	Used to contain red cells and logger if added								
Notes:	# If anticipated transport time exceeds the maximum transport time data logger must be placed in with the shipment.		Chilled ballast	2							

Lifeblood Shippers - Receipt and Use by External Institutions

R3	No. of Components per Shipper	Validated Transport Time [#]	Legend	Name	Quantity	Note																		
	1 to 14 red cell units	8hr 25min		Foil pouch	1	Used to line box and contains all other items																		
		Chilled ballast	4																					
		Extra chilled ballast (as per table) Place ballast outside the plastic liner bag	<table border="1"> <thead> <tr> <th>Number of red cell units</th> <th>Number of additional ballast packs</th> </tr> </thead> <tbody> <tr> <td>12 to 14</td> <td>0</td> </tr> <tr> <td>11</td> <td>1</td> </tr> <tr> <td>9 or 10</td> <td>2</td> </tr> <tr> <td>7 or 8</td> <td>3</td> </tr> <tr> <td>5 or 6</td> <td>4</td> </tr> <tr> <td>3 or 4</td> <td>5</td> </tr> <tr> <td>2</td> <td>6</td> </tr> <tr> <td>1</td> <td>7</td> </tr> <tr> <td>1 to 3 paediatric units</td> <td>7</td> </tr> </tbody> </table>		Number of red cell units	Number of additional ballast packs	12 to 14	0	11	1	9 or 10	2	7 or 8	3	5 or 6	4	3 or 4	5	2	6	1	7	1 to 3 paediatric units	7
	Number of red cell units	Number of additional ballast packs																						
	12 to 14	0																						
11	1																							
9 or 10	2																							
7 or 8	3																							
5 or 6	4																							
3 or 4	5																							
2	6																							
1	7																							
1 to 3 paediatric units	7																							
	Plastic liner bag	1	Used to contain red cells and logger if added																					
	Chilled ballast	4																						
Notes:	# If anticipated transport time exceeds the maximum transport time a data logger must be placed in with the shipment.																							

Lifeblood Shippers - Receipt and Use by External Institutions

R4	No. of Components per Shipper	Validated Transport Time [#]	Legend	Name	Quantity	Note																		
	Up to 10 red cell units	16hr 18min		Foil pouch	1	Used to line box and contains all other items																		
		Frozen -19°C	2																					
		Chilled ballast	2																					
		Cardboard divider	1	Used to stop contact with ballast																				
		Extra chilled ballast (as per table) Place ballast outside the plastic liner bag	<table border="1"> <thead> <tr> <th>Number of Red Cell Units</th> <th>Number of Ballast Packs</th> </tr> </thead> <tbody> <tr><td>10</td><td>0</td></tr> <tr><td>9</td><td>1</td></tr> <tr><td>7 or 8</td><td>2</td></tr> <tr><td>5 or 6</td><td>3</td></tr> <tr><td>3 or 4</td><td>4</td></tr> <tr><td>2</td><td>5</td></tr> <tr><td>1</td><td>6</td></tr> <tr><td>1 – 3 paediatric units</td><td>6</td></tr> <tr><td>1 – 3 Rh</td><td>6</td></tr> </tbody> </table>		Number of Red Cell Units	Number of Ballast Packs	10	0	9	1	7 or 8	2	5 or 6	3	3 or 4	4	2	5	1	6	1 – 3 paediatric units	6	1 – 3 Rh	6
	Number of Red Cell Units	Number of Ballast Packs																						
	10	0																						
	9	1																						
	7 or 8	2																						
	5 or 6	3																						
	3 or 4	4																						
2	5																							
1	6																							
1 – 3 paediatric units	6																							
1 – 3 Rh	6																							
	Plastic liner bag	1	Used to contain red cells and logger if added																					
	Cardboard divider	1	Used to stop contact with ballast																					
	Chilled ballast	2																						
	Frozen -19°C	2																						
Notes:	# If anticipated transport time exceeds the maximum transport time a data logger must be placed in with the shipment.																							

Lifeblood Shippers - Receipt and Use by External Institutions

4. Platelets

Platelet components should be transported using either **Packing Configuration P1 or P2**, the choice of which will depend on the anticipated transport time. If the transport time is estimated to approach or exceed the validated time, shown below, or if product is transported by independent couriers e.g. bus or aircraft, a data logger must be included.






Configuration P1 and P2 are designed to maintain the components within a temperature range of 20°C to 24 °C as recommended in the Council of Europe "Guide to the Preparation, Use and Quality Assurance of Blood Components" 14th edition.

Table 3 - Packing configurations and the validated transport time.








Packing configuration	No. of components per shipper	Validated transport time [#]
P1	1 to 7 pooled platelets 1-8 apheresis platelets 1 to 16 paediatric platelets	4hr 43min
P2	1 to 6 pooled platelets 1-8 apheresis platelets 1 to 16 paediatric platelets	7hr 38min
[#] If anticipated transport time exceeds the maximum transport time a data logger must be placed in the shipper between the components.		

Both configurations require a specific number of units to maintain temperature. Where fewer platelet units are to be packed, the specified number of **room temperature (+20 °C to +24 °C) ballast packs** must be added. The correct number of ballast packs for each configuration is shown in the packing diagrams below. The ballast packs should be placed outside the plastic liner to prevent contamination of the platelet units, should the ballast packs leak.

Lifeblood Shippers - Receipt and Use by External Institutions

P1	No. of Components per Shipper	Validated Transport Time[#]	Legend	Name	Quantity	Note																		
	1 to 7 pooled platelets 1 to 8 apheresis platelets 1 to 16 paediatric platelets	4 hr 43min		Room temperature ballasts 20-24°C	2																			
				Extra room temperature ballast (as per table) Place ballast outside the plastic liner bag	<table border="1"> <thead> <tr> <th colspan="2">Number of platelet units</th> <th>Number of additional ballast packs</th> </tr> </thead> <tbody> <tr> <td>4 - 8 (apheresis) 4 - 7 (pooled)</td> <td>Mixed total of 7</td> <td>0</td> </tr> <tr> <td>3</td> <td></td> <td>1</td> </tr> <tr> <td>2</td> <td></td> <td>2</td> </tr> <tr> <td>1</td> <td></td> <td>3</td> </tr> <tr> <td>1 to 16 paediatric</td> <td></td> <td>4</td> </tr> </tbody> </table>		Number of platelet units		Number of additional ballast packs	4 - 8 (apheresis) 4 - 7 (pooled)	Mixed total of 7	0	3		1	2		2	1		3	1 to 16 paediatric		4
			Number of platelet units		Number of additional ballast packs																			
4 - 8 (apheresis) 4 - 7 (pooled)	Mixed total of 7	0																						
3		1																						
2		2																						
1		3																						
1 to 16 paediatric		4																						
Notes:	# If anticipated transport time exceeds the maximum transport time a data logger must be placed in with the shipment.			Plastic liner bag	1	Used to contain platelets and logger if added																		
				Room temperature ballasts 20-24°C	2																			

Lifeblood Shippers - Receipt and Use by External Institutions

P2	No. of Components per Shipper	Validated Transport Time [#]	Legend	Name	Quantity	Note												
	1 to 6 pooled platelets 1 to 8 apheresis platelets 1 to 16 paediatric platelets	7hr 38min		Room temperature ballast 20-24°C	4													
 <p data-bbox="380 1109 683 1244">View of the long side of the shipper, placing two ballasts on each side.</p>				Foil pouch	1	Used to line box and contains all other items												
				Room temperature ballast 20-24°C	6	Place 2 along each of the long sides and 1 at each short side of the box												
				Extra room temperature ballast (as per the table) Place ballast outside the plastic liner bag but inside the foil pouch	<table border="1"> <thead> <tr> <th colspan="2">Number of platelet units</th> <th>Number of additional ballast packs</th> </tr> </thead> <tbody> <tr> <td>2 - 8 (apheresis) 2 - 6 (pooled)</td> <td>Mixed total of 7</td> <td>0</td> </tr> <tr> <td>1</td> <td></td> <td>1</td> </tr> <tr> <td>1 to 16 paediatric</td> <td></td> <td>2</td> </tr> </tbody> </table>		Number of platelet units		Number of additional ballast packs	2 - 8 (apheresis) 2 - 6 (pooled)	Mixed total of 7	0	1		1	1 to 16 paediatric		2
			Number of platelet units		Number of additional ballast packs													
			2 - 8 (apheresis) 2 - 6 (pooled)	Mixed total of 7	0													
1		1																
1 to 16 paediatric		2																
	Plastic liner bag	1	Used to contain platelets and logger if added															
	Room temperature ballast 20-24°C	4																

Notes: [#] If anticipated transport time exceeds the maximum transport time a data logger must be placed in with the shipment.

Configuration P2 requires 14 room temperature ballasts: 4 placed on the bottom, 2 on each of the longer sides, 1 on each of the shorted sides and another 4 ballast placed on the top.






5. Frozen components (plasma and cryoprecipitate)

The F2 Shipper configuration is used by Lifeblood for routine transport of frozen components. There are specific steps that must be taken when transporting shippers containing dry ice by road. Please refer to Appendix B for further information.





The F1 shipper configuration is available to use in the event that no dry ice is available, however it is not recommended for routine transport. If the use of this configuration is required, a data logger must be included in the shipment.

Configuration F1 and F2 are designed to maintain components at -25 °C as recommended in the Council of Europe "Guide to the Preparation, Use and Quality Assurance of Blood Components" 14th edition.

Lifeblood Shippers - Receipt and Use by External Institutions

	Number of components per polystyrene lined shipper	Validated Transport Time	Legend	Name	Quantity	Note											
F1	1 to 5 WB cFFP or APH cFFP	45min		Foil pouch	1	Used to line box and contains all other items											
	5 – 25 cryoprecipitate (no packaging)	*Refer to note below															
				Frozen - 40°C ballast	2	Place inside foil pouch – face down											
						<table border="1"> <thead> <tr> <th data-bbox="1464 603 1832 647">No. of Components per Shipper</th> <th data-bbox="1832 603 1968 647">Min</th> <th data-bbox="1968 603 2110 647">Max</th> </tr> </thead> <tbody> <tr> <td data-bbox="1464 647 1832 687">FFP whole blood / apheresis</td> <td data-bbox="1832 647 1968 687">1</td> <td data-bbox="1968 647 2110 687">5</td> </tr> <tr> <td data-bbox="1464 687 1832 727">Cryoprecipitate</td> <td data-bbox="1832 687 1968 727">5</td> <td data-bbox="1968 687 2110 727">25</td> </tr> </tbody> </table>			No. of Components per Shipper	Min	Max	FFP whole blood / apheresis	1	5	Cryoprecipitate	5	25
						No. of Components per Shipper	Min	Max									
			FFP whole blood / apheresis	1	5												
Cryoprecipitate	5	25															
	Plastic liner bag	1	Used to contain frozen components and logger if added														
	Frozen - 40°C ballast	2	Place inside foil pouch														
Notes:	*Use of configuration F1 is not recommended for the routine transport of cryoprecipitate. If use of this configuration is required a data logger must be included in all shipments.																

Lifeblood Shippers - Receipt and Use by External Institutions

F2	No. of Components per Shipper	Validated Transport Time [#]	Legend	Name	Quantity	Note																												
See table insert on right				Dry Ice	1.8 Kg	On top plastic liner pouch																												
				<table border="1"> <thead> <tr> <th>Number of Components</th> <th>Min</th> <th>Max</th> <th>Validated Transport Time[#]</th> </tr> </thead> <tbody> <tr> <td>Cryoprecipitate (no packaging)</td> <td>1</td> <td>25</td> <td>24hr</td> </tr> <tr> <td>Cryoprecipitate APH (no packaging)</td> <td>1</td> <td>8</td> <td>19hr:10min</td> </tr> <tr> <td>Cryodepleted APH 750ml (no packaging)</td> <td>1</td> <td>3</td> <td>19hr</td> </tr> <tr> <td>Cryoprecipitate (cardboard box)</td> <td>5</td> <td>25</td> <td>18hr:30min</td> </tr> <tr> <td>Cryoprecipitate APH (cardboard box)</td> <td>1</td> <td>6</td> <td>17hr:50min</td> </tr> <tr> <td>WB cFFP or APH cFFP (cardboard box & vacuum sealed)</td> <td>1</td> <td>6</td> <td>20hr:40min</td> </tr> </tbody> </table>			Number of Components	Min	Max	Validated Transport Time [#]	Cryoprecipitate (no packaging)	1	25	24hr	Cryoprecipitate APH (no packaging)	1	8	19hr:10min	Cryodepleted APH 750ml (no packaging)	1	3	19hr	Cryoprecipitate (cardboard box)	5	25	18hr:30min	Cryoprecipitate APH (cardboard box)	1	6	17hr:50min	WB cFFP or APH cFFP (cardboard box & vacuum sealed)	1	6	20hr:40min
Number of Components	Min	Max	Validated Transport Time [#]																															
Cryoprecipitate (no packaging)	1	25	24hr																															
Cryoprecipitate APH (no packaging)	1	8	19hr:10min																															
Cryodepleted APH 750ml (no packaging)	1	3	19hr																															
Cryoprecipitate (cardboard box)	5	25	18hr:30min																															
Cryoprecipitate APH (cardboard box)	1	6	17hr:50min																															
WB cFFP or APH cFFP (cardboard box & vacuum sealed)	1	6	20hr:40min																															
Notes: [#] If anticipated transport time exceeds the maximum transport time a data logger must be placed in with the shipment.				Plastic liner bag	1	Used to contain frozen components and logger if added																												
				Dry Ice	1.8 Kg	At bottom of shipper																												

Lifeblood Shippers - Receipt and Use by External Institutions

Appendix A: Details of suppliers of shipper consumables

Item	Company	Phone
Pink plastic bags (shipper liners)	Valpak	(02) 9984 0777
Shipper carton (outer cardboard box)	OJI Fibre Solutions	(02) 9724 8408
Shipper divider	OJI Fibre Solutions	(02) 9724 8408
Foil bags	Qualtape Australia	(03) 9729 8401
Polystyrene (Foam) inserts	R Max Rigid Cellular Plastics	(02) 9609 6088
Ballast pack with bubble wrap on one side	Sancell Pty Ltd	(03) 8796 5555
Temprecord data logger	Temprecord (NZ)	+64 9 274 9825
TempTale data logger	Sensitech	(03) 9686 5622

Appendix B: Transport of blood products using dry ice

The use of couriers must be taken into consideration as the use of dry ice in transporting items by vehicle can be hazardous to the driver, who may be unaware of the potential risk of exposure to CO₂ gas. Drivers need to be aware of the hazards and what to do in an emergency.

Below are two references which may be helpful for formulating policies and procedures.

- www.iata.org
- Technical Manual American Association of Blood Banks (AABB), 14th edition, pg 659

Lifeblood Shippers - Receipt and Use by External Institutions

Definitions

Term/abbreviation	Definition
Ballast pack	A gel pack that is preconditioned to maintain temperature for the transportation of blood products. Extra ballast packs (in excess of the base number of ballast required for a configuration) are added to provide temperature stability during transportation.
Blood products	The result of a procedure in which blood is withdrawn from a donor, separated into individual components e.g. platelets and/or retained as individual products for transfusion/fractionation.

Referenced external documents

N/A

Referenced internal documents

Document number	Document title
SOP-00326	Issuing and Packing Components and Products

Change history

Version number	Effective date	Reference	Summary of change
-	-	-	For previous change history contact National Document Control.
5	24/05/2021	CCR-20-000805	TempTale data logger supplier details added
6	12/07/2021	-	Minor grammatical changes to align with corresponding PPM configurations
		CCR-20-000781	Update validated transport times for R1, R2, R3, P1, P2 Update unit and ballast numbers for R1, R3. Remove R4 configuration.
7	16/07/2021	CCR-21-000600	Reinstate R4 configuration Updated R1 shipper to have VTT of 6 hours and not to be used for air travel
8	Refer to footer	CCR-21-000600	R4 – add '1' to number of red cell units.

Electronic signature

Author	NANCY KAY (nxkay)
Approver(s)	WAYNE ZHAO (wzhao)
Final Approver	MELANIE MADDOCK (mmaddock)
Owner	MELANIE MADDOCK (mmaddock)